

CERT Triage Pack Bandaging Project

Each triage pack needs (in gallon zip lock bag)

- 6-7 4" roller bandages
- 4-5 30"x30" triangular bandages


Washing Instructions

1. Wash sheets
2. Wash again in very hot water with bleach (according to package directions) with no detergent or fabric softener
3. Dry on hot in dryer
4. Fold and put in clean plastic bag

Cutting Instructions

1. Make sure the surface you are cutting on is clean!
2. Cut queen or king size sheets first. Cut queen or king size sheets into 4 triangular bandages and 4" wide roller bandages (see below). You should get approximately 10 bandages.
3. Cut smaller sheets. For smaller sheets, cut 4" roller bandages and/or triangular bandages, depending on need.
4. Roll bandages (and fold triangular bandages) and place in gallon zip lock bags.

Give to your Neighborhood Emergency Preparedness Director to put into neighborhood supplies.



EXAMPLE - Sheet sizes will dictate size of triangular bandages

CERT Triage Pack Bandaging Project

Each triage pack needs (in gallon zip lock bag)

- 6-7 4" roller bandages
- 4-5 30"x30" triangular bandages


Washing Instructions

5. Wash sheets
6. Wash again in very hot water with bleach (according to package directions) with no detergent or fabric softener
7. Dry on hot in dryer
8. Fold and put in clean plastic bag

Cutting Instructions

5. Make sure the surface you are cutting on is clean!
6. Cut queen or king size sheets first. Cut queen or king size sheets into 4 triangular bandages and 4" wide roller bandages (see below). You should get approximately 10 bandages.
7. Cut smaller sheets. For smaller sheets, cut 4" roller bandages and/or triangular bandages, depending on need.
8. Roll bandages (and fold triangular bandages) and place in gallon zip lock bags.

Give to Neighborhood Emergency Preparedness Director to put into neighborhood supplies.



EXAMPLE - Sheet sizes will dictate size of triangular bandages